



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Laboratory Metallurgical Microscope

## Product Image

---



## Description

---

<b>Body:</b>	Monocular inclined body.
<b>Stage:</b>	125x125mm
<b>Focussing:</b>	Separate fine and coarse focussing using diagonally cut rack and pinion arrangement.
<b>Illumination:</b>	6V 20Watt halogen illumination, operating on 220/110V, 50 Hz
<b>Eyepiece:</b>	Huygenian Eyepieces 5X, 10X with micrometer and 15X.
<b>Objective:</b>	M 4X/5X, M 10X, M 45X/40X

**Filter:** Day light blue filter.

**Box:** Supplied in thermocol packing.

Model	Catalog No.
Laboratory Metallurgical Microscope	JM006000

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.