



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Lactometer

Product Image





Description

Material:-Borosilicate Glass

A type Lactometers, Zeal Type

B type Lactometers, Quevene Type

C type Lactometers, with inbuilt Thermometer

A type: Lactometers, Zeal Type

Catalog No.	Specifications
--------------------	-----------------------

JG28101	Lactometer 0-40° Graduation at 84°F Zeal Type
JG28102	Lactometer 0-60° Graduation at 84°F Zeal Type
JG28103	Lactometer 0-40° Graduation at 84°F Zeal Type with mercury polish
JG28104	Lactometer 0-60° Graduation at 84°F Zeal Type with mercury polish

B type: Lactometers, Quevene Type

Catalog No.	Specifications
--------------------	-----------------------

JG28201	Lactometer Quevene type, 0-40° Graduation at 60°F
JG28202	Lactometer Quevene type, 15-40° Graduation at 60°F
JG28203	Lactometer Quevene type, 20-40° Graduation at 60°F
JG28204	Lactometer Quevene type, 15-35° Graduation at 60°F
JG28205	Lactometer Quevene type, 20-35° Graduation at 80°F

C type: Lactometers, with inbuilt Thermometer

Catalog No.	Specifications
--------------------	-----------------------

JG28301	Lactometer 0-40°, with inbuilt Thermometer 0-110°C, Graduation at 60°F
JG28302	Lactometer 20-35°, with inbuilt Thermometer 0-40°, Graduation at 68°F

JG28303

Lactometer 20-40°, with inbuilt Thermometer
0-40°, Graduation at 68°F

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.