




Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

MCT Twin Rack

Product Image



Description

Material: Polypropylene

Provides this twin Purpose rack for both 0.5ml & 1.5ml. Micro centrifge tubes. Each side has 108 holes for 0.5ml & 1.5ml MCTs respectively. The transparent top provides an excellent view as well as a protection cover for the tubes stores. Each face of the base has a unique alpha-numeric index for easy indentification of the tube.

Catalog No.	Particulars
JW05201	Rack for 0.5ml and 1.5ml MCTs

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.