



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Masses (Weights), Hooked, Double sided, Recessed with Sliding Hook, Brass

## Product Image



## Description

---

Polished brass masses (weights) with hook at both ends. The bottom is recessed and hooks have sliding mechanism. This enables suspension from both ends and enables masses be placed on flat surface. In wooden block.

Catalog No.	Model
-------------	-------

JW76800	Capacity 2100g (1x10g, 1x20g, 1x20g, 1x50g, 1x100g, 1x200g, 1x200g, 1x500g & 1x1000g)
---------	---

**Note:** Individual masses (weights) and custom-made sets available on request.

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.