



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Masses (Weights), Slotted in Rack

Product Image



Description

Precision polished brass masses (weights). These slotted masses can be placed on a rack for both horizontal and vertical storage.

**Catalog
No.**

Capacity(in gm)

JW77301

Capacity 1100g (1x10g,
2x20g,1x50g,1x100g, 2x200g, &
1x500g)

JW77302 Capacity 1110g (1x1g,
2x2g,1x5g,1x10g,
2x20g,1x50g,1x100g, 2x200g &
1x500g)

Note: Individual masses (weights) and custom-made sets available on request.

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.