



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Masses (Weights) Set, Slotted, Brass

Product Image



Description

Precision brass masses (weights) and hangers, Hanger made of brass wire with hook for suspension. Specially designed slot of the masses to prevent the masses from slipping off the hangers.

Catalog No.	Capacity(in gm)
JW77401	One hanger of 5g and nine slotted masses of 5g each, total 50g
JW77402	One hanger of 10g and nine slotted masses of 10g each, total 100g
JW77403	One hanger of 20g and nine slotted masses of 20g each, total 200g
JW77404	One hanger of 50g and nine slotted masses of 50g each, total 500g
JW77405	One hanger of 100g and nine slotted masses of 100g each, total 1000g
JW77406	One hanger of 200g and nine slotted masses of 200g each, total 2000g
JW77407	One hanger of 20g, three slotted masses at 20g each, one slotted mass of 10g, and two masses of 5g each, total 100g
JW77408	One hanger of 50g. nine slotted masses of 20g each, one slotted mass of 10g, and two slotted masses of 5g each, total 250g

Note: Individual masses (weights), hangers and custom-made sets available on request.

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.