



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Masses (Weights) Set, Steel, Hooked, Painted

Product Image



Description

Complete set of painted hooked masses (weights) made of Steel for general purpose weighing. In wooden block.

Catalog No.	Capacity(in gm)
JW77201	Capacity 100g (1x50g, 2x20g,1x10g)
JW77202	Capacity 200g (1x100g, 1x50g, 2x20g, 1x10g)
JW77203	Capacity 500g (1x200g, 2x100g, 1x50g, 2x20g, 1x10g)
JW77204	Capacity 1000g (1x500g,1x200g, 2x100g,1x50g, 2x20g, 1x10g)
JW77205	Capacity 2000g (1x1000g,1x500g,1x200g, 2x100g,1x50g, 2x20g, 1x10g)

Note: Individual masses (weights) and custom-made sets available on request.

Note: Also available in Black CED (Cathode Electro Deposition) finish

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.