



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Masses set (Weights), Knobbed, Brass

## Product Image

---



## Description

---

Precision masses (weights) made of polished brass in knobbed design. In wooden block.

<b>Catalog No.</b>	<b>Capacity( in gm)</b>
JW76901	Capacity 100g (1x50g, 1x20g, 2x10g,1x5g, 2x2g,1x1g)
JW76902	Capacity 200g (1x100g, 1x50g, 1x20g, 2x10g,1x5g, 2x2g,1x1g)
JW76903	Capacity 500g (1x200g, 2x100g, 1x50g, 1x20g, 2x10g,1x5g, 2x2g,1x1g)
JW76904	Capacity 1000g (1x500g,1x200g,2x100g,1x50g,1x20g, 2x10g,1x5g, 2x2g,1x1g)
JW76905	Capacity 2000g (1x1000g,1x500g,1x200g, 2x100g,1x50g,1x20g, 2x10g,1x5g, 2x2g,1x1g)

**Note:** Following masses (weights) sets are also available in Brass with Chrome Plating

**Note:** Individual masses (weights) and custom-made sets available on request.

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.