



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Megenta Box

## Product Image

---



## Description

---

**Material:** Polycarbonate/Polypropylene

Magenta box is very often used in tissue culture and agriculture research experiments. Main body of this box is made of tough and transparent polycarbonate which provides excellent view of the contents inside while the cover is made of polypropylene. This autoclavable box is square in shape giving itself a space saving feature.

**Catalog No.**

JW12301

**Item**

Magenta Box

## Disclaimer

---

The Product details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.