



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Micro Centrifuge, 13000 R.P.M.

## Product Image



## Description

---

Model : T-13

### Specifications:

- High speed Micro Centrifuge Max. R.P.M. 13000.
- Heavy duty plastic injection moulded body.
- Stainless Steel protection bowl.
- Fitted with well balance universal motor having long life with high grade carbons.
- 2 lines of 16 character LCD display of r.p.m., r.c.f, set time and run time.
- Digital timer which can be set from 0-99 minutes in interval of 1 minute.
- Programmable Speed Regulator from 1000 to 13000 in variation of 100 rpm.
- Last test parameter recall and displays automatically on LCD.

### Supplied completed with following accessories:

- Rotor Head to accommodate 24 eppendorf/Micro Centrifuge Tubes with cover of 1.5ml /2ml (24 x 1.5/2ml).
- Reduction adopters for 0.2/0.5ml eppendorf (Micro Centrifuge Tubes).
- Supply is complete with head Puller, Spanner, Instruction manual, Cord & plug to work on 220 volts 50Hz A.C.

-

Dimension of Machine: L = 290mm, W = 270mm, H = 240mm
Weight of Machine: 9.5kg. approx.
Dimension of box: L-355mm x W-355 x H-340mm
Weight of Box :11 kg. approx.

-

Micro Centrifuge, 13000 R.P.M.

**Catalog No. 100764**

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.