



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Micro Tip Box

Product Image



Description

Material: Polypropylene

These unique box are moulded in polypropylene and are therefore autoclavable. The box itself is clear providing an excellent view of the tips inside while the table like structure is provided inside to hold the tips. The tip holding space are uniformly placed to facilitate the use of multi channel pipettes. The cover is joined to the base through ball catchers hinges and provides excellent protection to the tips.

| Catalog No. | Particulars |
|--------------------|-----------------------------------|
| JW07701 | For 96 Micro Tips of 2-200 ul |
| JW07702 | For 100 Micro Tips of 200-1000 ul |

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.