




**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Micro Tip Box

## Product Image



## Description

---

**Material : PP**

Micro tip boxes molded in high quality polypropylene can be autoclavable. The unique feature of these tip boxes is the robust hinge design which prevents the top from being easily disengaged from the bottom of the box and consistent wall thickness for repetitive use in single and multi channel pipette. Transparent top cover provides excellent view of the micro tips inside the box.

Catalog No.	Particulars
101445-A	For 96 Tips of 2-10µl
101445-B	For 96 Tips of 2-20µl
101445-C	For 100 Tips of 200-1000µl

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.