



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Micropipette, Regular, Fixed Volume

Product Image



Description

Micropipette, Regular is a compact device with separate tip ejector. These micropipettes enjoys excellent reputation with medical researchers and professional people. It comes with superior features to achieve highest efficiency. The unit is user friendly and strain free. Each pipette is calibrated to specifications conforming to DIN 12650 Standards. These pipettes are partially autoclavable.

Features :

Colour coded for easy identification
New ergonomic design
It reduced probability of repetitive strain injury
Separate tip ejector
Clear digital window on the hand grip

Particulars

Catalog No.

Micropipette, Regular, Fixed Volume 2 ul	JA90520
Micropipette, Regular, Fixed Volume 5 ul	JA90521
Micropipette, Regular, Fixed Volume 10 ul	JA90522
Micropipette, Regular, Fixed Volume 20 ul	JA90523
Micropipette, Regular, Fixed Volume 25 ul	JA90524
Micropipette, Regular, Fixed Volume 50 ul	JA90525
Micropipette, Regular, Fixed Volume 100 ul	JA90526
Micropipette, Regular, Fixed Volume 200 ul	JA90527
Micropipette, Regular, Fixed Volume 250 ul	JA90528
Micropipette, Regular, Fixed Volume 500 ul	JA90529
Micropipette, Regular, Fixed Volume 1000 ul	JA90530
Micropipette, Regular, Fixed Volume 2000 ul	JA90531
Micropipette, Regular, Fixed Volume 5000 ul	JA90532

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.