



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Microprocessor Haemoglobin Meter

Product Image



Description

Catalog No: 100102

Microprocessor Haemoglobin Meter - Model 113

TECHNICAL SPECIFICATIONS

✓ DISPLAY	3 DIGIT RED LED DISPLAY.
✓ DISPLAY MODES	% T, O.D. AND 3 CONCENTRATION MODES C1,C2,C3
✓ ZERO ADJUSTMENT	AUTOMATIC (WITHOUT KNOB ROTATION)
✓ ACCURACY	± 0.02 O.D.
✓ STABILITY	± 0.02 O.D. PER HOUR.
✓ LIGHT SOURCE	6.8 VOLT , 0.3 AMP TUNGSTEN LAMP.
✓ AMPLIFIER	MICRO CONTROLLER BASED.
✓ MEMORY	IN BUILT MEMEORY TO RETAIN LAST "K" FACTOR FOR CONCENTRATION EVEN AFTER POWER LOSS.
✓ SAMPLE VOLUME	1 ML. MINIMUM.
✓ POWER SUPPLY	230 V \pm 10 %, 50 Hz.

ACCESSORIES

✓ Glass test tubes	5 Nos.
✓ Instruments cover	1 Nos.
✓ Operating Manual	1 Nos.
✓ Spare Lamp	1 Nos. (Fitted inside)

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.