



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Microprocessor pH Meter, Superior (5 Point Calibration)

Product Image



Description

Catalog No : 100093

Microprocessor pH Meter, Superior (5 Point Calibration) - Model 50

TECHNICAL SPECIFICATIONS

PARAMETER	pH	mV	Temp^oC
✓ RANGE	0.00 to 14.00 pH	0 to 1999.9 mV	0 to 100 ^o C
✓ RESOLUTION	0.01pH	0.1 mV	0.1 ^o C
✓ ACCURACY	±0.01 pH ±1 digit	±0.1mV	±0.1 ^o C
✓ TEMP. COMPENSATION	-5 to 110 °C		
✓ CALIBRATION	Auto/Manual		
✓ AUTO BUFFER	4.00, 7.00 & 9.20 pH		
✓ READ OUT	128 x 64 dots Graphical LCD Display		
✓ STORAGE	Up to 1000 Samples		
✓ PC ATTACHMENT	Through USB		
✓ KEY BOARD	5 Keys Soft Touch Membrane Type		
✓ POWER	230±10% AC		
✓ WEIGHT	0.8 Kg. (Approx.)		

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.