



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Mini-PCR

Product Image



Description

Features:

- It is small-sized and easy to program with an intuitive user interface.
- The lid adopts the high temperature resistance material and application to various types of test tubes.
- Memory function incase of power-down.
- PC and K160 can be separately controlled through the PC operating software.
- It is beneficial for students to understand with animation presentation capabilities of the PC operating software.

Parameter

Capacity	25×0.2mL(A); 9×0.5mL(B); 16×0.2mL+9×0.5mL(C)
Temp range	0℃~99.9℃(Rt≤30℃)
Max heating rate	≥2℃/s
Max cooling rate	≥2℃/s
Uniformity	≤±0.3℃ (Constant 10s)
Accuracy	≤±0.2℃
Temp control	Block
Stored program No.	3
Max. No. of Cycle	99
Display	12864LCD
Communication	serial port
Size	160mm(L)x140mm(W)x120mm(H)
Weight	2.2kg

-
Mini-PCR

Catalog No. 101050

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.