



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Monocular Research Microscope, Basic

## Product Image



## Description

---

<b>Body:</b>	45 degrees inclined monocular head rotatable through 360 degrees. Quadruple nose piece with accurate centering and positive click stops.
<b>Stage:</b>	Mechanical stage of 120x140mm. X-Y movement with fully graduated vernier scale for manipulation of specimen slide.
<b>Focussing:</b>	Separate fine and coarse focussing using diagonally cut rack and pinion arrangement.
<b>Illumination:</b>	6V-20W halogen lamp, operating on 220/110 Volts.
<b>Condenser:</b>	Two lens Abbe condenser N.A. 1.25 with filter support, movable by rack & pinion arrangement.
<b>Eyepiece:</b>	Wide Field 10X and either of Huygenian 5X or 15X.
<b>Objective:</b>	Achromatic JIS Standard 36mm 4X, 10X, 40X SL and 100X SL oil immersion
<b>Box:</b>	Supplied in a plywood box with essential accessories.

Model	Catalog No.
Monocular Research Microscope, Basic with eyepiece Wide Field 10X, Huygenian 5X	JM001201
Monocular Research Microscope, Basic with eyepiece Wide Field 10X, Huygenian 15X	JM001202

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.