




Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Orsat Gas Analysis Apparatus 3 Test

Product Image



Description

Material:- Borosilicate Glass

**Catalog
No.**

Capacity (ml)

JG19501	Orsat Gas Analysis Apparatus 3 Tests: For the determination of CO ₂ , CO and O ₂ particularly in flue and furnace gases, comprising of levelling bottle, 100ml Gas burette with outer jacket, three absorption pipettes, manifold with Glass key stopcocks in wooden case with sliding doors.
JG19502	With manifold having Glass key stopcock similar to "Orsat Gas Analysis Apparatus 3 Test" but with four absorption pipettes.
JG19503	With manifold having Glass key stop cock similar to "Orsat Gas Analysis Apparatus 3 Test" in design but with four absorption pipettes and a palladium /platinum asbestos tube. This is for the analysis of CO ₂ , CO, O, H ₂ and nitrogen (By Difference).
JG19504	With Five absorption pipettes for determination of CO ₂ , CO, H ₂ , O ₂ . Illuminants methane, ethane and nitrogen by difference. Comprising of four absorption pipettes, two of them plain type fitted with Glass tubes, two bobbling type pipettes with three way

Note: All the above models are available with PTFE (Teflon) Key Stopcock

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.