



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Orthopedic OT Table

## Product Image



# Description

---

Model 504

-

## (Hydraulic C ARM Orthopedic OT Table)

- | Radio translucent top in four sections for use of C arm image intensifier.
- | Compatible for orthopedic attachments.
- | Heavy duty smooth gear mechanism to avail all operating positions.
- | Durable heavy duty leak proof hydraulic pump with four pillars.
- | Eccentric top and hydraulic chamber adjusted at one side of base for maximum use of C ARM image intensifier.
- | Changeable oil valve for up & down movement of hydraulic pump.
- | Inter changeable head and leg section.
- | One radio translucent hand board and one hand rest with S.S. top for transfusion.
- | Most useful for orthopedic surgery.
- | Specially designed for intra operative X ray & C arm image intensifier.
- | X ray translucent removable pelvis rest with special arrangement for orthopedic operation.
- | Padded orthopedic knee crutches provides support to legs of patient.
- | Two traction device sliding support with ball joints allowing flxion & rotations of foot with traction shoes.
- | Provided with hand traction attachment.

## STANDARD ACCESSORIES

- | Anesthetic screen.
- | Shoulder supports with pads.
- | Arm boards with pads.
- | Side support with pads
- | Black rubberized mattress.

-

Orthopedic OT Table                      **Catalog No. 100230**

# Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.