



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

P.R.P. Test Centrifuge , Type 2

Product Image



Description

Model : 8M

Specifications:

- Compact P.R.P Test Centrifuge Machine, Microprocessor based, digital.
- Body made of rust proof ABS Plastic injection moulding with metallic protection bowl.
- Lid lock for safe drive.
- 2 lines of 16 character LCD display of r.p.m., for set and run time.
- Two Programs are feeded in software of machine as per recommended by for P.R.P. Tests.
 - (i) First Program 4 minutes, 3400 r.p.m.

(ii) Second Program 2 minutes, 3500 r.p.m. OR
we can programme speed from 500 to 5000 r.p.m. in variation of 100 r.p.m. as per Doctors requirement .

- Last test parameter recall and displays automatically on LCD.
- Sample Holding Capacity: 4 P.R.P. Syringes of 20 ml.
- Dr. P.R.P. Kit: Having 20 ml syringe can be used with our machine.

Technical Parameters:

Motor	Universal Motor (with carbon brushes) 1/8 H.P.
Power	AC 220V 50Hz $\pm 10\%$
Dimension of Machine	L = 415mm, W = 360mm, H = 305mm
Dimension of outer box	510 x 510 x 390mm
Net Weight	21 kg. Approx. Gross Weight with box = 23.5 Kg.
Current	220 volts, 50Hz A.C.

P.R.P. Test Centrifuge , Type 2

Catalog No. 100755

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.