



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

PCR Plate

Product Image



Description

PCR Plates are moulded from high quality virgin polypropylene in 96 well format. The thin wall V-bottom design provides rapid and efficient thermal transfer to each sample. These PCR Strips are compatible with most thermal cyclers.

Catalog No.	Capacity
JW07881	56 well non skirted, elevated well
JW07882	96 well skirted

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.