



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

PCR Tube Rack

Product Image



Description

Material: Polycarbonate/Polypropylene This rack is for 96 PCR tubes of 0.2ml, as the name suggests, is used for keeping PCR tubes. Top & base of this rack are made of tough & transparent polycarbonate while the middle portion is made of polypropylene. The tube holding spaces are uniformly places so as to facilitate the use of multipipettes . A no. of racks can be stacked together to save shelf space.

Catalog No.	Particulars
JW05301	Rack for 96 PCR tubesof 0.2ml

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.