



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

PCR Tube Strip

Product Image



Description

PCR Strips are manufactured from high quality virgin polypropylene and have thin uniform wall for efficient thermal transfer. These PCR Strips are designed for optimal PCR performance.

Catalog No.	Capacity (ml)
JW07877	0.2 ml strip of 8
JW07878	0.2 ml strip of 12

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.