



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# PCR Tubes

## Product Image



## Description

---

Flat cap PCR tubes are manufactured from high quality virgin polypropylene and have thin uniform wall for efficient thermal transfer. These PCR tubes are designed for optimal performance.

Catalog No.	Capacity (ml)
JW07875	0.2
JW07876	0.2 Dnase/Rnase Free Certified

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.