




Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

PIPETTE PUMP

Product Image



Description

25 ml Rubber Pipette Pump A simple, manual pipette/dropper/syringe with a transparent glass body for easy visibility. It features a straight or curved tip, with a rubber bulb at the opposite end. The pipette is graduated, with capacities of 0.75 ml, 0.50 ml, 0.25 ml, and a total capacity of 1 ml up to 25 ml.

Catalog no	Particulars
------------	-------------

103216	PIPETTE PUMP
--------	--------------

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.