



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# PP Trays for Lab Furniture

## Product Image


---



## Description

---

We offers 100% moulded Polypropylene Trays which are made from highly chemical resistant PP material which are available in various sizes and are used in various applications such as Fume Hood, Reagent Rack and Chemical Storage Unit.



Length (mm)	Width (mm)	Height (mm)
755	745	45
975	335	45
895	335	45
795	335	45
625	335	45
405	335	45

PP Trays for Lab Furnitures

Catalog No. 101173

## Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.