



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Penta Head Multi Viewing Microscope, Research

Product Image



Description

Multi-viewing Microscope allows people observation at the same time. With the features of infinite optical system, effective illumination, LED pointer and images. Coherence, it is widely used in clinic, scientific research and teaching demonstration.

Body:	Consist of compensation free 1 piece trinocular head and 4 pieces binocular heads, inclined at 30°, rotatable at 360°, interpupillary distance 48-75mm with inverse quintuple nosepiece.
Stage:	Double layers mechanical stage 185x142mm, with moving range 75x55mm.
Focussing:	Standard co-axial coarse and fine focussing with fine division upto 0.002mm.
Illumination:	Kohler external illumination, halogen lamp 12V/50W.
Eyepiece:	Extra wide field eyepiece EW10X/20.
Objective:	Infinite plan achromatic objective 4X,10X, 40X and 100X.
Condenser:	Swing condenser.
Pointer:	Green LED pointer, brightness adjustable.
Box:	Supplied in thermocol packing. Accessories:

Accessories:

1. Spare halogen lamp
2. Immersion oil
3. Dust cover
4. User manual
5. Blue ground glass filter

Optional Accessories

Pointer: Dual color pointer, brightness adjustable.

Photo Attachment

Video Adapter:

- C Mount 1X
- C Mount 0.5X

Model

Catalog No.

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.