



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Petri Dish Box, Aluminium/Stainless Steel

Product Image



Description

Petri dish box is seamless, die pressed box with a removable steel rack to hold 100mm Ø

petri plates.

Protects Petri dishes during and after heat sterilization.

Polished Stainless Steel, including deep telescoping cover and removable rack.

Available in Stainless steel and Aluminium.

Petri Dish Box, Aluminium, 4 inches x 10 inches

Catalog No. 100201-A

Petri Dish Box, Stainless Steel, 4 inches x 10 inches

Catalog No. 100201-B

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.