



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Petri Dishes Disposable - Gamma Sterile

Product Image



Description

Material: Polystyrene

Petri dishes are made of polystyrene. Its dimensional accuracy of 90mm ensures compatibility with most automated systems. Its optical clarity ensures a clear view on microscope. These triple vented dishes have an elevated lid that allows for good gas exchange & encourages evaporation. These petri dishes are sterilized by radiation.

-

Catalog No.

Particulars

101469 - A

Petri Dishes 90 mm (individually packed) - Gamma Sterile

101469 - B

Petri Dishes (Sleeve of 20) - Gamma Sterile

-

Note: Also available Petri Dishes Disposable ETO Sterile

Note: Also available Petri Dishes Disposable Aseptic

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.