



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Physiograph Single Channel

Product Image



Description

Digital Student Physiograph with Single Channel Console with Time & Event Channel and stimulator for Human & Animal experiments.

Main Console with Time & Event Channel and stimulator. High Resolution coloured TFT display Sweep speeds provided are 0.5,1,2,5,10,20,25,30,50 div/sec. The data channel has recording width of 8 div.(80mm approx.)

Couplers:-

- i) Strain Gage
- ii) Isotonic
- iii) Pulse Respiration
- iv) Temperature
- v) EKG (Clinical) with Electrode, 5 pin Junction box and Jelly.
- vi) Bio Potential (with Electrodes, 3 pin Junction box, Paste and Electrode for Action Potential).

Transducers:-

- i) Pressure for strain gage Coupler
- ii) Volume for strain gage Coupler
- iii) Pulse for Pulse Respiration Coupler
- iv) Temperature for Temperature Coupler
- v) Muscle Activity/Force for strain gage Coupler
- iv) Respiration Belt for Pulse Respiration Coupler
- vii) Isotonic Fine Movement for isotonic Coupler
- viii) Respiration (Thermister Type) for Pulse Respiration Coupler

It is also available in 2 channel, 3 channel and 4 channel.

-

Physiograph Single Channel

Catalog No. 100662

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.