



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Pipette Box, Aluminium/Stainless Steel

Product Image



Description

Pipette Box is Designed to secure and protect pipette during sterilization processes.

Made for Stainless Steel or Aluminium construction to ensures high corrosion resistance. Ideal for preventing breakage of costly pipettes and glass tubes during sterilization either by autoclaving or dry heat.

Multiple box styles accommodate a variety of pipette types.

Pipette Box is avaialable in various sizes:

Pipette Box, Stainless Steel, 2.5 inches x 16 inches	Catalog No. 100202-A
Pipette Box, Stainless Steel, 3 inches x 18 inches	Catalog No. 100202-B
Pipette Box, Aluminium, 2.5 inches x16 inches	Catalog No. 100202-C

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.