



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Pipette Serological

Product Image



Description

Pipette Serological, are manufactured exclusively with high grade polystyrene(GPPS) which ensures excellent clarity. The filter prevents overfilling and suitable with all brands of

pipetting apparatus. Extremely Bold. Black graduations with clear scale ensure accuracy.

Catalog No.	Capacity (ml)
JW13341	1
JW13342	2
JW13343	5
JW13344	10
JW13345	25
JW13346	50

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.