



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Pipette Stand (Vertical)

Product Image



Description

Material: Polypropylene

This is a stable non corrosive detachable stand which can hold 28 pipettes with maximum dia 14 mm. The heavy base measures 225 mm in dia and a polypropylene coated square rod measuring 375 mm in length is screwed into the base. The top disc that can be fixed on the rod has holes for holding pipettes while the lower one has pin holes where the pipette tips rest and liquids drain down to the base. These discs can be adjusted on the rod accordingly to the length of the pipette. The arrows provided on the disc should be kept parallel to ensure that corresponding holes are being centered.

| Catalog No. | Particulars |
|-------------|-----------------|
| JW09501 | For 28 Pipettes |

Disclaimer

The Product details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.