



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Pipette Tips & Fliter Tips

Product Image



Description

Material: Polypropylene

Catalog No.	Types	Capacity(µl)
JW13311	Clear, Non Sterilized	0.1-10
JW13312	Clear, Sterilized	0.1-10
JW13313	Clear, Sterilized + Certified Dnase/Rnase Free	0.1-10
JW13314	Clear, with Filter Non Sterilized	0.1-1
JW13315	Clear, with Filter Sterilized	0.1-10
JW13316	Clear, with Filter Sterilized + Certified Dnase/Rnase Free	0.1-10
JW13317	Yellow, Non Sterilized	10-200
JW13318	Yellow, Sterilized	10-200
JW13319	Clear, Sterilized + Certified Dnase/Rnase Free	10-200
JW13320	Clear, with Filter Sterilized + Certified Dnase/Rnase Free	10-200
JW13321	Blue, Non Sterilized	100-1000
JW13322	Blue, Sterilized	100-1000
JW13323	Clear, Sterilized + Certified Dnase/Rnase Free	100-1000
JW13324	Clear, with Filter Non Sterilized	100-1000
JW13325	Clear, with Filter Sterilized	100-1000
JW13326	Clear, with Filter Sterilized + Certified Dnase/Rnase Free	100-1000

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.