



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Plasme Freezer

Product Image



Description

Outer body is made up of thick mild sheet and inner chamber is of stainless steel. For uniform and safe storage of plasma and related blood components at uniform preset temperature range - 30°C. All vertical models of plasma freezers are fitted with 7 days circular chart recorder fitted Digital Temp. indicator - cum Controller supplied with three

stainless steel storing shelves. To work on 220 volts.

Also Available:

- Digital temperature controller /PID/TOUCH SCREEN LCD controller with data logging and computer interfacing options on Demand.
- Complete unit in 304 grade stainless steel on request.

Number of Bags:

-130 BAGS

-240 BAGS

Plasme Freezer

Catalog No. 100396

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.