



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Portable Haemoglobin Meter

## Product Image



## Description

Catalog No : 100103

Portable Haemoglobin Meter - Model 113P

Portable and high performance digital hemoglobin meter is requires fully automated process. This provides easy, fast & very accurate monitoring of haemoglobin. In this the adjustment or calibration is not needed and ensures to provide results within 15 seconds.

## TECHNICAL SPECIFICATIONS

- ✓ 10 Testing Strips
- ✓ Measuring range: 5-25 gm/dl
- ✓ Memory: 1000 test results
- ✓ Auto calibration

## ACCESSORIES

- ✓ HB Meter
- ✓ Canister of Test Strips
- ✓ Code Chip
- ✓ Control chip
- ✓ Capillary Trasfer Tube/Dropper
- ✓ Lancing Device
- ✓ Sterile Lancets
- ✓ AAA Batteries

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.