



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Portable Haemoglobin Meter

Product Image



Description

Catalog No : 100103

Portable Haemoglobin Meter - Model 113P

Portable and high performance digital hemoglobin meter is requires fully automated process. This provides easy, fast & very accurate monitoring of haemoglobin. In this the adjustment or calibration is not needed and ensures to provide results within 15 seconds.

TECHNICAL SPECIFICATIONS

- ✓ 10 Testing Strips
- ✓ Measuring range: 5-25 gm/dl
- ✓ Memory: 1000 test results
- ✓ Auto calibration

ACCESSORIES

- ✓ HB Meter
- ✓ Canister of Test Strips
- ✓ Code Chip
- ✓ Control chip
- ✓ Capillary Trasfer Tube/Dropper
- ✓ Lancing Device
- ✓ Sterile Lancets
- ✓ AAA Batteries

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.