



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Precision Rotary Microtome, Spencer Type

## Product Image

---



## Description

---

Most modern and dependable sectioning equipment specially designed for precise sectioning tissues upto 1 micron thickness. This model is a product of meticulous quality and embodies all the features of a standard sectioning instrument to give perfect and accurate results.

- Base is massive and rigid to provide stability.
- Protective cover is hinged to the base and swings back to completely expose the interior parts for easier cleaning and lubrication.
- Feed range : 1 to 50 microns in steps of 1 micron, precise feed setting with latest cam drive system.
- Automatic feed release, safety device and wheel lock.
- Feed Excursion : 35 mm.
- Opening of object clamp 30x15mm.
- Overall size : 40 x 25 x 23 cms. approx
- Supplied with complete accessories three object holder, one honing plate, one knife and one oil can for lubrication & Wooden cabinet.

### Particulars

### Catalog No.

Precision Rotary Microtome, Spencer  
Type

JM00450

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.