



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Precision Rotary Microtome, Spencer Type

Product Image



Description

Most modern and dependable sectioning equipment specially designed for precise sectioning tissues upto 1 micron thickness. This model is a product of meticulous quality and embodies all the features of a standard sectioning instrument to give perfect and accurate results.

- Base is massive and rigid to provide stability.
- Protective cover is hinged to the base and swings back to completely expose the interior parts for easier cleaning and lubrication.
- Feed range : 1 to 50 microns in steps of 1 micron, precise feed setting with latest cam drive system.
- Automatic feed release, safety device and wheel lock.
- Feed Excursion : 35 mm.
- Opening of object clamp 30x15mm.
- Overall size : 40 x 25 x 23 cms. approx
- Supplied with complete accessories three object holder, one honing plate, one knife and one oil can for lubrication & Wooden cabinet.

Particulars

Catalog No.

Precision Rotary Microtome, Spencer
Type

JM00450

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.