



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Projection Microscope with LCD Screen

Product Image



Description

Body:	Well designed body with sturdy base with the screen fitted on the tube of the microscope, rotatable in all direction.
Focussing:	Separate fine and coarse focussing using diagonally cut rack and pinion arrangement.
Illumination:	3W LED super white light illumination, operating on 220/110 Volts.
Box:	Supplied in thermocol packing.
Camera Specification:	
Image Sensor:	1/3.2" CMOS
Valid Pixel:	1600x1200 (2MP)
Pixel Size:	2.8 μ x2.8 μ
Digital Output:	24bit (color)
Image Format:	Frame rate:- 1600x1200 7.5f/s:800x600 30f/s
White Balance:	Manual/Auto white Balance
Output:	AV output

TFT LCD TV/Monitor

1. LCD Size: 9" TFT LCD Screen.
2. Eight languages menu changeable & blue background display.
3. All function controlled by remote control.
4. Multi System, PAL/NTSC/SECAM
5. AV Input.
6. DC 12V
7. USB Support LCD TFT Screen.

Model	Catalog No.
Projection Microscope with LCD Screen	JM002600

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.