




**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# RED ALLIGATOR CLIP TEST LEADS

## Product Image



## Description

Test leads with a length ranging from 30 to 50 cm, a width of at least 2 mm, and red color. They feature dual alligator clips made of copper material at both ends, offering flexibility. The wires are insulated with high-quality PVC for durability.

Catalog no	Particulars
103643	RED ALLIGATOR CLIP TEST LEADS

## Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.