



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Rack For Centrifuge Tube (Folding) - Space Saver

Product Image



Description

Material: Polypropylene

These space saver folding centrifuge tube racks are molded in PP. Have raised tube sockets which guarantees firm hold of the tubes at their exacts position. The alpha numeric grid reference on the top tier ensures easy tube location & sample positioning. They are supplied assembled and folded, can be ready for use in few seconds with a simple action and can be folded again after use for reducing storage area. These racks are autoclavable & can with stand temperature ranging between - 80°C to + 121°C. These racks are available in individual for 15ml to 50 ml tubes and in combination for both tubes.

Catalog No.

Particulars

101473 -A

Space Saver - 15ml- 45 Place

101473 -B

Space Saver - 50ml- 18 Place

101473 -C

Space Saver 15 & 50ml- 21 Place & 12 Place resp.

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.