



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com


**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Rack For Micro Centrifuge (Folding)

## Product Image

---



## Description

---

**Material:** Polypropylene

These two-tier versatile racks, made of polypropylene, are autoclavable and can be used dry, in water bath or in freezers. These racks are provided with alpha numeric index for easy sample identification. Racks are stackable even when loaded and have purpose designed label spaces on each side.

**Catalog No.**

**Particulars**

JW09101

100 Places for 1.5 ml MCTs

## Disclaimer

---

The Product details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.