



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Rack for 15/50 ml Centrifuge Tubes

Product Image



Description

This rack is most suitable where conical centrifuge tubes are in job. It has alphanumeric inscription for easy identification of test and can be easily assembled or disassembled. Withstands autoclaving, freezing.

Catalog No.

JW08621

Capacity (ml)

15 ml/50 places

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.