



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Rack for Micro Centrifuge Tubes

Product Image



Description

Material: Polycarbonate

Moulded in polycarbonate these autoclavable racks are made for holding 1.5 ml micro centrifuge tubes. These racks have a unique stacking features which prevent movement in any direction. The holes are numbered and lettered for easy identification.

Catalog No.	Particulars
JW08601	24 Tubes
JW08602	48 Tubes

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.