



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Rack for Micro Centrifuge Tubes

## Product Image

---



## Description

---

**Material:** Polycarbonate

Moulded in polycarbonate these autoclavable racks are made for holding 1.5 ml micro centrifuge tubes. These racks have a unique stacking features which prevent movement in any direction. The holes are numbered and lettered for easy identification.

<b>Catalog No.</b>	<b>Particulars</b>
JW08601	24 Tubes
JW08602	48 Tubes

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.