



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Racks For Petri Dishes

Product Image



Description

Material: Polymethyl Methacrylate/Polycarbonate

These racks are clear acrylic construction with poly carbonate posts 69101 has a capacity to hold 60 petri dish of 90 mm size while 69102 one can hold 56 petri dish of 60 mm. These racks are very useful during inoculation, incubation, and storage. Visibility from all sides enables the culture to be seen and checked during incubation.

Catalog No.	Particulars
JW08101	Racks of 90 mm Petri Dish
JW08102	Racks of 60 mm Petri Dish

Disclaimer

The Product details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.