



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Racks For Scintillation Vial

Product Image



Description

Material: Polypropylene

These racks for scintillation vials, moulded in polypropylene are autoclavable and corrosion resistant. There is an individual rack for both 8 ml and 20 ml vials. These rugged and stable racks can hold scintillation vials for convenient storage handling and sample preparation.

Catalog No.	Particulars
JW06901	Rack for 90 vials of 8 ml cap
JW06902	Racks for 50 vials of 20 ml cap

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.