



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Reagent Bottle (Narrow Mouth), Plastic

## Product Image

---



## Description

---

**Material :** Polypropylene These narrow Mouth Reagent Bottles are available in Polypropylene. The PP Bottles are rigid, translucent and have excellent conatct clarity. These bottles can be autoclaves filled as well as empty. The Polypropylene Cap has a built seal ring to make these bottles leakproof.

Catalog No.	Capacity (ml)
JW00801	4
JW00802	8
JW00803	15
JW00804	30
JW00805	60
JW00806	125
JW00807	250
JW00808	500
JW00809	1000

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.