



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Reagent Bottle (Wide Mouth), Plastic

Product Image



Description

Material : Polypropylene These Wide Mouth Reagent Bottles are available in Polypropylene. The PP bottles are rigid, translucent and have excellent contact clarity. These can be autoclaved filled as well as empty. The Polepropylene cap has a built in seal ring to make these bottles leak proof.

Catalog No.	Capacity (ml)
JW00901	30
JW00902	60
JW00903	125
JW00904	250
JW00905	500
JW00906	1000

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.